

Anti-Nuak1 Rabbit Polyclonal Antibody For Wb, Elisa - 060285

Item Number: CM0018934

Introduction

Rabbit polyclonal antibody targeting human NUA1/ARK5 kinase, validated for western blotting and ELISA applications. Detects both human and rat proteins; ideal for investigating cell adhesion, senescence, and cancer progression pathways.

[Learn More](#)

Parameter	Value
Product Name	NUAK1 Rabbit pAb
Remarks/Alias	ARK5, NUA1
Species	Human
GeneID (Human)	9891
GeneID	9891
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 362-661 of human NUA1 (NP_055655.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	O60285
Protein Weight	74kDa
Shipping	Ice bag