

Anti-Nuak2 Polyclonal Antibody For If/Icc - Q9H093

Item Number: CM0023960

Introduction

Rabbit polyclonal antibody against human NUA2 (SNARK), validated for immunofluorescence/immunocytochemistry. Used for research on apoptosis, cell migration, Hippo signaling, and neural tube development.

[Learn More](#)

Parameter	Value
Product Name	NUAK2 Rabbit pAb
Remarks/Alias	NUAK2; SNARK; NUA family kinase 2
Species	human
GenelD (Human)	81788
GenelD	81788
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC
Cross Reactivity	Human
SWISS	Q9H093
Protein Weight	69kDa
Shipping	Ice bag