


Anti-Nub1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y5A7

Item Number: CM0018792



Introduction

NUB1 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA in Human, Mouse, Rat. Targets NEDD8 ultimate buster 1, a down-regulator of the NEDD8 conjugation system, suitable for proteasome degradation research.

[Learn More](#)

Parameter	Value
Product Name	NUB1 Rabbit pAb
Remarks/Alias	BS4; NUB1L; NYREN18; NUB1
Species	Human
GeneID (Human)	51667
GeneID	51667
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human IFT140 (NP_001350458.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5A7
Protein Weight	71kDa
Shipping	Ice bag