

Anti-Nubp1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P53384

Item Number: CM0020059

Introduction

NUBP1 Rabbit Polyclonal Antibody, validated for WB, IHC-P, ELISA. Targets cytosolic Fe-S cluster assembly factor NUBP1, cross-reacts with Human, Mouse, Rat. Ideal for studying iron-sulfur protein biogenesis.

[Learn More](#)

Parameter	Value
Product Name	NUBP1 Rabbit pAb
Remarks/Alias	NBP; NBP1; CIAO5; NBP35; NUBP1
Species	Human
GeneID (Human)	4682
GeneID	4682
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 201-320 of human NUBP1 (NP_002475.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P53384
Protein Weight	35kDa
Shipping	Ice bag