

Anti-Nucb2 Polyclonal Antibody For Wb, Elisa - P80303

Item Number: CM0014253

Introduction

Rabbit polyclonal antibody targeting human NUCB2, validated for WB and ELISA. Suitable for calcium homeostasis, energy metabolism, and G protein signaling research. Cross-reacts with human and rat.

[Learn More](#)

Parameter	Value
Product Name	NUCB2 Rabbit pAb
Remarks/Alias	NEFA; HEL-S-109; NUCB2
Species	Human
GeneID (Human)	4925
GeneID	4925
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 321-420 of human NUCB2 (NP_005004.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P80303
Protein Weight	50kDa
Shipping	Ice bag