

Anti-Nucks1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9H1E3

Item Number: CM0022402

Introduction

High-quality NUCKS1 Rabbit Polyclonal Antibody validated for Western blot (WB) and ELISA, showing cross-reactivity with human, mouse, and rat. Ideal for studying DNA repair via homologous recombination and chromatin dynamics.

[Learn More](#)

Parameter	Value
Product Name	NUCKS1 Rabbit pAb
Remarks/Alias	JC7; NUCKS; NUCKS1
Species	Human
GeneID (human)	64710
GeneID	64710
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human NUCKS1 (NP_073568.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H1E3
Protein Weight	27kDa
Shipping	Ice bag