


Anti-Nuclear Matrix Protein P84 (Thoc1) Monoclonal Antibody For Wb, Ip, Elisa - Q96Fv9

Item Number: CM0011268



Introduction

Rabbit monoclonal antibody targeting Nuclear Matrix Protein p84 (THOC1) with SWISS accession Q96FV9. Suitable for western blot, immunoprecipitation, and ELISA in human, mouse, rat samples. Investigate TREX complex and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	Nuclear Matrix Protein p84 (THOC1) Rabbit mAb
Remarks/Alias	P84; HPR1; P84N5; DFNA86; Nuclear Matrix Protein p84 (THOC1)
Species	Human
Gene ID (Human)	9984
Gene ID	9984
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 350-450 of human Nuclear Matrix Protein p84 (THOC1) (NP_005122.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96FV9
Protein Weight	76kDa
Shipping	Ice bag