


Anti-Nuclear Matrix Protein P84 (Thoc1) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q96Fv9

Item Number: CM0027912



Introduction

Rabbit polyclonal antibody against human THOC1/p84 (76kDa). Suitable for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse, rat. Explores mRNA export, apoptosis, cancer.

[Learn More](#)

Parameter	Value
Product Name	Nuclear Matrix Protein p84 (THOC1) Rabbit pAb
Remarks/Alias	P84; HPR1; P84N5; DFNA86; Nuclear Matrix Protein p84 (THOC1)
Species	Human
GeneID (human)	9984
GeneID	9984
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 428-657 of human Nuclear Matrix Protein p84 (NP_005122.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96FV9
Protein Weight	76kDa
Shipping	Ice bag