

Anti-Nucleobindin-1 (Nucb1) Rabbit Polyclonal Antibody For Wb And Elisa - Q02818

Item Number: CM0021541



Introduction

Rabbit polyclonal antibody targeting human nucleobindin-1 (NUCB1/CALNUC), validated for Western blot and ELISA. Suitable for investigating Golgi calcium-binding protein functions and G protein signaling.

[Learn More](#)

Parameter	Specification
Product Name	NUCB1 Rabbit pAb
Remarks/Alias	NUC; CALNUC; NUCB1
Species	Human
GeneID (Human)	4924
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 302-461 of human NUCB1 (NP_006175.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q02818
Protein Weight	54kDa
Shipping	Ice bag