

Anti-Nucleolin Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P19338

Item Number: CM0002390



Introduction

Rabbit monoclonal antibody targeting human Nucleolin (C23), applicable for Western blot, IHC-P, IF/ICC, and ELISA. Suitable for chromatin and ribosome biogenesis research.

[Learn More](#)

Parameter	Value
Product Name	Nucleolin Rabbit mAb
Remarks/Alias	C23; Nsr1; Nucleolin
Species	Human
Gene ID	4691
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P19338
Protein Weight	77kDa
Shipping	Ice bag