

Anti-Nucleolin (Ncl) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P19338

Item Number: CM0025732



Introduction

High-quality anti-NCL rabbit polyclonal antibody validated for WB, IHC-P, ELISA. Recognizes human and mouse nucleolin, a major nucleolar protein. Ideal for chromatin and ribosome research.

[Learn More](#)

Parameter	Value
Product Name	NCL Rabbit pAb
Remarks/Alias	C23; Nsr1; NCL
Species	Human
Gene ID (Human)	4691
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-466 of human NCL (NP_005372.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P19338
Protein Weight	77kDa
Shipping	Ice bag