

Anti-Nucleophosmin Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P06748

Item Number: CM0008538



Introduction

Anti-Nucleophosmin monoclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. For ribosome biogenesis, centrosome duplication, and tumor suppressor regulation research.

[Learn More](#)

Parameter	Value
Product Name	Nucleophosmin Rabbit mAb
Remarks/Alias	B23; NPM
Species	Human
GeneID (Human)	4869
GeneID	4869
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 9-122 of human NPM1 (NP_001238905.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06748
Protein Weight	33kDa
Shipping	Ice bag