


Anti-Nudc Polyclonal Antibody For Wb, Elisa - Q9Y266

Item Number: CM0027519



Introduction

Rabbit polyclonal antibody against human NUDC (Nuclear migration protein nudC), validated for WB and ELISA. Suitable for detecting NUDC in human, mouse, and rat samples. Recombinant immunogen, 38 kDa.

[Learn More](#)

Parameter	Value
Product Name	NUDC Rabbit pAb
Remarks/Alias	HNUDC; MNUDC; NPD011; NUDC
Species	Human
GeneID (Human)	10726
GeneID	10726
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-331 of human NUDC (NP_006591.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y266
Protein Weight	38kDa
Shipping	Ice bag