


Anti-Nudc Rabbit Monoclonal Antibody For Wb And Elisa - Q9Y266

Item Number: CM0009737



Introduction

High-quality rabbit monoclonal antibody targeting human NUDC, suitable for Western blot and ELISA applications. Cross-reacts with mouse and rat NUDC. Ideal for studying mitosis, cytokinesis, and neurogenesis.

[Learn More](#)

Parameter	Value
Product Name	NUDC Rabbit mAb
Remarks/Alias	HNUDC; MNUDC; NPD011; NUDC
Species	Human
GeneID (Human)	10726
GeneID	10726
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 232-331 of human NUDC (Q9Y266).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y266
Protein Weight	38kDa
Shipping	Ice bag