

Anti-Nudel/Ndel1 Monoclonal Antibody For Wb, Elisa - Q9Gzm8

Item Number: CM0010254



Introduction

Anti-NUDEL/NDEL1 rabbit monoclonal antibody validated for Western blot and ELISA; reacts with human, mouse, rat. Suitable for research in neurodevelopment, cytoskeletal organization, and mitosis.

[Learn More](#)

Parameter	Value
Product Name	NUDEL/NDEL1 Rabbit mAb
Remarks/Alias	EOPA; NDE2; NUDEL; MITAP1; NDE1L1; NUDEL/NDEL1
Species	Human
Gene ID (Human)	81565
Gene ID	81565
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human NUDEL/NDEL1 (Q9GZM8).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9GZM8
Protein Weight	38kDa
Shipping	Ice bag