

Anti-Nudt15 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Nv35

Item Number: CM0013311



Introduction

Rabbit polyclonal antibody against human NUDT15, suitable for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for studying nucleotide metabolism, thiopurine drug response, and DNA synthesis.

[Learn More](#)

Parameter	Value
Product Name	NUDT15 Rabbit pAb
Remarks / Alias	MTH2; NUDT15D; NUDT15
Species	Human
GeneID (Human)	55270
GeneID	55270
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-164 of human NUDT15 (NP_060753.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9NV35
Protein Weight	19kDa
Shipping	Ice bag