

Anti-Nudt16 Polyclonal Antibody For Wb, Elisa - Q96De0

Item Number: CM0021023



Introduction

Rabbit polyclonal antibody against NUDT16, a U8 snoRNA-decapping enzyme and nucleotide pool sanitizer. Validated for WB and ELISA in human and mouse. Ideal for RNA metabolism and nucleic acid quality control studies.

[Learn More](#)

Parameter	Value
Product Name	NUDT16 Rabbit pAb
Remarks/Alias	NUDT16
Species	Human
Gene ID	131870
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 60-130 of human NUDT16 (NP_689608.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96DE0
Protein Weight	21kDa
Shipping	Ice bag