

Anti-Nudt18 Polyclonal Antibody For Wb, Elisa - Q6Zvk8

Item Number: CM0028240



Introduction

Rabbit polyclonal antibody targeting human NUDT18 (MTH3), validated for Western blot and ELISA. Suitable for DNA/RNA damage and nucleotide metabolism research.

[Learn More](#)

Parameter	Value
Product Name	NUDT18 Rabbit pAb
Remarks/Alias	MTH3; NUDT18
Species	Human
GeneID (Human)	79873
GeneID	79873
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-323 of human NUDT18 (NP_079091.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q6ZVK8
Protein Weight	36kDa
Shipping	Ice bag