

Anti-Nudt2 Polyclonal Antibody For Wb, Elisa - P50583

Item Number: CM0026982



Introduction

NUDT2 rabbit polyclonal antibody validated for WB and ELISA. Recognizes human, mouse, and rat NUDT2. Ideal for studying Ap4A hydrolase function, RNA decapping, and related neurological disorders.

[Learn More](#)

Parameter	Value
Product Name	NUDT2 Rabbit pAb
Remarks/Alias	APAH1; IDDPN; NUDT2
Species	Human
GeneID (Human)	318
GeneID	318
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-147 of human NUDT2 (NP_001152.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50583
Protein Weight	17kDa
Shipping	Ice bag