

Anti-Nudt21 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O43809

Item Number: CM0029014



Introduction

NUDT21 rabbit polyclonal antibody for Western blot, IHC, and ELISA. Detects human, mouse, rat NUDT21. Ideal for studying pre-mRNA cleavage and polyadenylation, and alternative polyadenylation research.

[Learn More](#)

Parameter	Value
Product Name	NUDT21 Rabbit pAb
Remarks/Alias	CPSF5; CFIM25; NUDT21
Species	Human
GeneID (Human)	11051
GeneID	11051
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-227 of human NUDT21 (NP_008937.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43809
Protein Weight	26kDa
Shipping	Ice bag