


Anti-Nudt5 Monoclonal Antibody For Wb, If/Icc, Elisa - Q9Ukk9

Item Number: CM0002278



Introduction

Rabbit monoclonal antibody targeting human NUDT5, suitable for WB, IF/ICC, and ELISA applications. Cross-reacts with human, mouse, and rat. NUDT5 is an ADP-sugar pyrophosphatase involved in nuclear ATP synthesis.

[Learn More](#)

Parameter	Value
Product Name	NUDT5 Rabbit mAb
Remarks/Alias	YSA1; YSA1H; YSAH1; hNUDT5; NUDT5
Species	Human
GeneID (Human)	11164
GeneID	11164
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human NUDT5 (Q9UKK9).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UKK9
Protein Weight	24kDa
Shipping	Ice bag