

Anti-Nudt6 Polyclonal Antibody For Wb, If/Icc, Elisa - P53370

Item Number: CM0018311



Introduction

NUDT6 Rabbit pAb is a polyclonal antibody targeting human NUDT6, validated for WB, IF/ICC, and ELISA. NUDT6 hydrolyzes NADH and FAD and may regulate cell proliferation, making it valuable for nucleotide metabolism and mitochondrial studies.

[Learn More](#)

Parameter	Value
Product Name	NUDT6 Rabbit pAb
Remarks/Alias	GFG1; GFG-1; ASFGF2; FGF-AS; FGF2AS; NUDT6
Species	Human
GeneID (Human)	11162
GeneID	11162
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 121-316 of human NUDT6 (NP_009014.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P53370
Protein Weight	36kDa
Shipping	Ice bag