


Anti-Nuf2 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Bzd4

Item Number: CM0003012



Introduction

High-quality rabbit monoclonal anti-NUF2 antibody, validated for WB and ELISA applications. Reacts with human NUF2, a kinetochore protein essential for chromosome segregation; predicted molecular weight 54kDa. Ideal for research in cell cycle and mitosis.

[Learn More](#)

Parameter	Value
Product Name	NUF2 Rabbit mAb
Remarks/Alias	CDCA1; CT106; NUF2R; NUF2
Species	Human
GeneID	83540
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human NUF2 (Q9BZD4).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9BZD4
Protein Weight	54kDa
Shipping	Ice bag