


Anti-Numa1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q14980

Item Number: CM0015396



Introduction

Rabbit polyclonal antibody to human NUMA1 (Nuclear Mitotic Apparatus Protein 1, NMP-22). Validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for studying mitotic spindle dynamics and nuclear matrix organization.

[Learn More](#)

Parameter	Value
Product Name	NUMA1 Rabbit pAb
Remarks/Alias	NUMA; NMP-22; NUMA1
Species	Human
Gene ID (Human)	4926
Gene ID	4926
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-309 of human NUMA1 (NP_006176.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14980
Protein Weight	238kDa
Shipping	Ice bag