


Anti-Numb Rabbit Monoclonal Antibody For Wb, Ip, Elisa - P49757

Item Number: CM0011421



Introduction

High-quality anti-NUMB rabbit monoclonal antibody validated for WB, IP, and ELISA in human and mouse samples. Targets protein numb homolog (P49757), a key regulator of endocytosis and neurogenesis.

[Learn More](#)

Parameter	Value
Product Name	NUMB Rabbit mAb
Remarks/Alias	S171; C14orf41; c14_5527; NUMB
Species	Human
GeneID (human)	8650
GeneID	8650
Immunogen	A synthetic peptide corresponding to aa 1-100 of human NUMB (P49757).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P49757
Protein Weight	71kDa
Shipping	Ice bag