

Anti-Numb Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P49757

Item Number: CM0027769



Introduction

Rabbit polyclonal antibody targeting human NUMB, suitable for WB, IHC-P, IF/ICC, and ELISA. Detects endogenous NUMB in human, mouse, and rat samples. Ideal for neurogenesis and endocytosis research.

[Learn More](#)

Parameter	Value
Product Name	NUMB Rabbit pAb
Remarks/Alias	S171; C14orf41; c14_5527; NUMB
Species	Human
Gene ID	8650
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 521-640 of human NUMB (NP_003735.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49757
Protein Weight	71kDa
Shipping	Ice bag