

Anti-Nup107 Polyclonal Antibody For Wb, Elisa - P57740

Item Number: CM0013201



Introduction

Rabbit polyclonal antibody against human NUP107, validated for Western blot (WB) and ELISA. Detects NUP107 in human, mouse, and rat samples. NUP107 is a nuclear pore complex protein essential for NPC assembly and nephrogenesis.

[Learn More](#)

Parameter	Value
Product Name	NUP107 Rabbit pAb
Remarks/Alias	ODG6; NUP84; NPHS11; GAMOS7; NUP107
Species	Human
Gene ID (human)	57122
Gene ID	57122
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 766-925 of human NUP107 (NP_065134.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P57740
Protein Weight	106kDa
Shipping	Ice bag