

Anti-Nup153 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P49790

Item Number: CM0006793



Introduction

CamelBio's NUP153 rabbit monoclonal antibody detects human, mouse, and rat nucleoporin 153 (154 kDa) for Western blotting, immunohistochemistry, immunofluorescence, and ELISA. Ideal for studying nuclear pore complex, nucleocytoplasmic transport, and viral host interactions.

[Learn More](#)

Parameter	Value
Product Name	NUP153 Rabbit mAb
Remarks/Alias	N153; HNUP153; NUP153
Species	Human
GeneID (Human)	9972
GeneID	9972
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human NUP153 (P49790)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49790
Protein Weight	154kDa
Shipping	Ice bag