

Anti-Nup155 Rabbit Monoclonal Antibody For Wb, Elisa - O75694

Item Number: CM0004416



Introduction

Anti-NUP155 rabbit monoclonal antibody validated for Western blot and ELISA. Recognizes human and mouse NUP155, a 155 kDa nuclear pore complex protein essential for nucleocytoplasmic transport. Ideal for cell biology and transport mechanism studies.

[Learn More](#)

Parameter	Value
Product Name	NUP155 Rabbit mAb
Remarks/Alias	N155; ATFB15; NUP155
Species	Human
GeneID (Human)	9631
GeneID	9631
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 901-1050 of human NUP155 (NP_705618.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O75694
Protein Weight	155kDa
Shipping	Ice bag