

Anti-Nup155 Rabbit Polyclonal Antibody For Wb, Elisa - 075694

Item Number: CM0025416



Introduction

Rabbit polyclonal antibody targeting human NUP155 (155 kDa nucleoporin), validated for Western blot (WB) and ELISA. Essential nuclear pore complex component for nucleocytoplasmic transport; ideal for IVD research and assay development.

[Learn More](#)

Parameter	Value
Product Name	NUP155 Rabbit pAb
Remarks/Alias	N155; ATFB15; NUP155
Species	Human
Gene ID (Human)	9631
Gene ID	9631
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1102-1391 of human NUP155 (NP_705618.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O75694
Protein Weight	155kDa
Shipping	Ice bag