

Anti-Nup160 Rabbit Polyclonal Antibody For Wb, Elisa - Q12769

Item Number: CM0030666



Introduction

NUP160 Rabbit Polyclonal Antibody validated for WB and ELISA with human reactivity. Targets nuclear pore complex protein Nup160, a 162kDa nucleoporin involved in mRNA transport. Ideal for NPC and nucleocytoplasmic transport research.

[Learn More](#)

Parameter	Value
Product Name	NUP160 Rabbit pAb
Remarks/Alias	NPHS19; NUP160
Species	Human
GeneID	23279
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-150 of human NUP160 (NP_056046.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q12769
Protein Weight	162kDa
Shipping	Ice bag