


Anti-Nup205 Rabbit Polyclonal Antibody For Wb And Elisa - Q92621

Item Number: CM0023641



Introduction

High-quality rabbit polyclonal antibody targeting human NUP205, validated for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for investigating nuclear pore complex function and nucleocytoplasmic transport mechanisms.

[Learn More](#)

Parameter	Value
Product Name	NUP205 Rabbit pAb
Remarks/Alias	NPHS13; C7orf14; NUP205
Species	Human
Gene ID (Human)	23165
Gene ID	23165
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 10-120 of human NUP205 (NP_055950.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92621
Protein Weight	228kDa
Shipping	Ice bag