


Anti-Nup210 (Gp210/Pom210) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Tem1

Item Number: CM0021245



Introduction

Rabbit polyclonal antibody targeting NUP210 (GP210, POM210), a 205 kDa nuclear pore membrane glycoprotein. Validated for WB, IF/ICC, and ELISA, with cross-reactivity to human and mouse samples. Suitable for nuclear pore complex assembly and transport research.

[Learn More](#)

Parameter	Value
Product Name	NUP210 Rabbit pAb
Remarks/Alias	GP210; POM210; NUP210
Species	Human
GeneID (human)	23225
GeneID	23225
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 680-880 of human NUP210 (NP_079199.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8TEM1
Protein Weight	205kDa
Shipping	Ice bag