


Anti-Nup214 Polyclonal Antibody For Wb, If/Icc, Elisa - P35658

Item Number: CM0030911



Introduction

Rabbit polyclonal antibody targeting human NUP214 (CAN protein). Validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for studying nucleocytoplasmic transport and nuclear pore complexes.

[Learn More](#)

Parameter	Value
Product Name	NUP214 Rabbit pAb
Remarks/Alias	CAN; CAIN; IIAE9; NUP214
Species	Human
GeneID (human)	8021
GeneID	8021
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1251-1350 of human NUP214 (NP_005076.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35658
Protein Weight	214kDa
Shipping	Ice bag