


Anti-Nup35 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Nfh5

Item Number: CM0013808



Introduction

Rabbit polyclonal antibody targeting human nucleoporin NUP35, validated for Western blot and ELISA. Reacts with human, mouse, and rat samples. Suitable for nuclear pore complex and transport research.

[Learn More](#)

Parameter	Value
Product Name	NUP35 Rabbit pAb
Remarks/Alias	MP44; NP44; MP-44; NUP53; NUP35
Species	Human
Gene ID (Human)	129401
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-326 of human NUP35 (NP_612142.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NFH5
Protein Weight	35kDa
Shipping	Ice bag