


Anti-Nup50 Rabbit Polyclonal Antibody - Q9Ukx7

Item Number: CM0027310



Introduction

Rabbit polyclonal antibody targeting human NUP50 (Q9UKX7) for WB, IHC-P, ELISA applications. Validated for human, mouse, and rat samples; suitable for nuclear pore complex and protein import studies.

[Learn More](#)

Parameter	Value
Product Name	NUP50 Rabbit pAb
Remarks/Alias	NPAP60; NPAP60L; NUP50
Species	Human
GeneID (Human)	10762
GeneID	10762
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UKX7
Protein Weight	50kDa
Shipping	Ice bag