


Anti-Nup54 Rabbit Polyclonal Antibody For Wb, Elisa - Q7Z3B4

Item Number: CM0015308



Introduction

Anti-NUP54 rabbit polyclonal antibody validated for WB and ELISA. Detects human, mouse, and rat NUP54 protein (~55 kDa). Ideal for nuclear pore complex and nucleocytoplasmic transport research.

[Learn More](#)

Parameter	Value
Product Name	NUP54 Rabbit pAb
Remarks/Alias	NUP54
Species	Human
Gene ID (Human)	53371
Gene ID	53371
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 100-400 of human NUP54 (NP_059122.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7Z3B4
Protein Weight	55kDa
Shipping	Ice bag