


Anti-Nup62 Monoclonal Antibody For Wb, Ip, Elisa - P37198

Item Number: CM0002200



Introduction

Anti-NUP62 rabbit monoclonal antibody validated for WB, IP, and ELISA; reacts with human, mouse, and rat NUP62. Targets nuclear pore complex protein involved in nucleocytoplasmic transport and mitotic progression. For research use only.

[Learn More](#)

Parameter	Value
Product Name	NUP62 Rabbit mAb
Remarks/Alias	p62, IBSN, SNDI, NUP62
Species	Human
GeneID (Human)	23636
GeneID	23636
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 423-522 of human NUP62 (P37198)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P37198
Protein Weight	53kDa
Shipping	Ice bag