

Anti-Nup85 Polyclonal Antibody For Wb, Elisa - Q9Bw27

Item Number: CM0013040



Introduction

Rabbit polyclonal antibody against human NUP85, a nuclear pore complex protein involved in NPC assembly, RNA export, mitosis, and monocyte chemotaxis. Validated for WB and ELISA; cross-reacts with human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	NUP85 Rabbit pAb
Remarks/Alias	Nup75; FROUNT; NPHS17; NUP85
Species	Human
GeneID (human)	79902
GeneID	79902
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 427-656 of human NUP85 (NP_079120.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BW27
Protein Weight	75kDa
Shipping	Ice bag