


Anti-Nup93 Rabbit Polyclonal Antibody For Wb, Elisa - Q8N1F7

Item Number: CM0027379



Introduction

NUP93 rabbit polyclonal antibody for Western blot (WB) and ELISA. Detects human, mouse, rat NUP93. Suitable for nuclear pore complex assembly and podocyte migration research.

[Learn More](#)

Parameter	Specification
Product Name	NUP93 Rabbit pAb
Remarks/Alias	NIC96; NUP93
Species	Human
GeneID (human)	9688
GeneID	9688
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 670-819 of human NUP93 (NP_055484.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N1F7
Protein Weight	93kDa
Shipping	Ice bag