

Anti-Nup98 Polyclonal Antibody For Wb, If/Icc, Elisa - P52948

Item Number: CM0014393



Introduction

This rabbit polyclonal antibody targets NUP98, a nuclear pore complex protein involved in nucleocytoplasmic transport. Validated for WB, IF/ICC, and ELISA applications, it recognizes human and mouse NUP98, making it suitable for NPC assembly and trafficking studies.

[Learn More](#)

Parameter	Value
Product Name	NUP98 Rabbit pAb
Remarks/Alias	ADIR2; NUP96; NUP196; Nup98-96; NUP98
Species	Human
GeneID (Human)	4928
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 621-920 of human NUP98 (NP_624357.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P52948
Protein Weight	198kDa
Shipping	Ice bag