

Anti-Nup98 Rabbit Monoclonal Antibody For Wb - P52948

Item Number: CM0007086



Introduction

NUP98 Rabbit mAb validated for WB and ELISA. Detects endogenous NUP98 across human, mouse, and monkey species. Essential tool for nuclear pore complex transport and assembly studies.

[Learn More](#)

Parameter	Value
Product Name	NUP98 Rabbit mAb
Remarks/Alias	ADIR2; NUP96; NUP196; Nup98-96; NUP98
Species	Human
Gene ID	4928
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 712-863 of human NUP98 (NP_001352055.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Monkey
SWISS	P52948
Protein Weight	198kDa
Shipping	Ice bag