


## Anti-Nur77 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P22736

Item Number: CM0024922



### Introduction

Rabbit polyclonal antibody targeting human NUR77 (NR4A1) validated for WB, IHC-P, ELISA. Reacts with human and mouse. Useful for studying nuclear receptor signaling, inflammation, and apoptosis. 64kDa protein, shipped on ice bags.

[Learn More](#)

| Parameter        | Specification                                                                                                              |
|------------------|----------------------------------------------------------------------------------------------------------------------------|
| Product Name     | NUR77 Rabbit pAb                                                                                                           |
| Remarks/Alias    | HMR; N10; TR3; NP10; GFRP1; NAK-1; NGFIB; NUR77                                                                            |
| Species          | Human                                                                                                                      |
| GeneID (Human)   | 3164                                                                                                                       |
| GeneID           | 3164                                                                                                                       |
| Immunogen        | Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 300-400 of human NUR77 (NP_002126.2). |
| Source           | Rabbit                                                                                                                     |
| Category         | Polyclonal Antibodies                                                                                                      |
| Application      | WB, IHC-P, ELISA                                                                                                           |
| Cross Reactivity | Human, Mouse                                                                                                               |
| SWISS            | P22736                                                                                                                     |
| Protein Weight   | 64kDa                                                                                                                      |
| Shipping         | Ice bag                                                                                                                    |