

Anti-Nusap1 Polyclonal Antibody For Wb And Elisa - Q9Bxs6

Item Number: CM0019201

Introduction

Rabbit polyclonal antibody targeting NUSAP1 (aa 1-61). Validated for Western blot and ELISA. Reacts with human and mouse NUSAP1. Ideal for mitosis and cell cycle research.

[Learn More](#)

Parameter	Value
Product Name	NUSAP1 Rabbit pAb
Remarks/Alias	LNP; ANKT; SAPL; BM037; NUSAP; Q0310; PRO0310p1; NUSAP1
Species	Human
Gene ID (Human)	51203
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-61 of human NUSAP1 (NP_060924.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BXS6
Protein Weight	49kDa
Shipping	Ice bag