

## Anti-Nutf2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P61970

Item Number: CM0013620



### Introduction

High-quality rabbit polyclonal antibody against human NUTF2 (NTF2/PP15), validated for WB, IF/ICC, and ELISA. Ideal for nucleocytoplasmic transport research. Cross-reacts with human and rat.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | NUTF2 Rabbit pAb                                                                                                 |
| Remarks/Alias    | NTF2; PP15; NTF-2; NUTF2                                                                                         |
| Species          | Human                                                                                                            |
| Gene ID (Human)  | 10204                                                                                                            |
| Gene ID          | 10204                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-127 of human NUTF2 (NP_005787.1) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, IF/ICC, ELISA                                                                                                |
| Cross Reactivity | Human, Rat                                                                                                       |
| SWISS            | P61970                                                                                                           |
| Protein Weight   | 14kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |