

Anti-Nxf1 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q9Ubu9

Item Number: CM0008727



Introduction

High-quality anti-NXF1 (TAP) rabbit monoclonal antibody validated for human WB, IP, and ELISA. Targets nuclear RNA export factor 1, a key mRNA export receptor. Predicted molecular weight 70 kDa.

[Learn More](#)

Parameter	Value
Product Name	NXF1 Rabbit mAb
Remarks/Alias	TAP; MEX67; NXF1
Species	Human
GeneID (Human)	10482
GeneID	10482
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-92 of human NXF1 (Q9UBU9)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	Q9UBU9
Protein Weight	70kDa
Shipping	Ice bag