

Anti-Nxf1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Ubu9

Item Number: CM0029590



Introduction

Rabbit polyclonal antibody targeting human NXF1 (TAP/MEX67), suitable for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for mRNA export and nuclear transport research.

[Learn More](#)

Parameter	Value
Product Name	NXF1 Rabbit pAb
Remarks/Alias	TAP; MEX67; NXF1
Species	Human
GeneID (Human)	10482
GeneID	10482
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 200-300 of human NXF1 (NP_006353.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBU9
Protein Weight	70kDa
Shipping	Ice bag