

Anti-Nxt2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Npj8

Item Number: CM0015571



Introduction

Rabbit polyclonal antibody targeting human NXT2 (NTF2-related export protein 2). Validated for WB, IF/ICC, and ELISA with cross-reactivity to mouse and rat. Essential tool for studying mRNA and protein nuclear export.

[Learn More](#)

Parameter	Value
Product Name	NXT2 Rabbit pAb
Remarks/Alias	P15-2; NXT2
Species	Human
GeneID (Human)	55916
GeneID	55916
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 110-180 of human NXT2 (NP_061168.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NPJ8
Protein Weight	16kDa
Shipping	Ice bag