

Anti-Ny-Eso-1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P78358

Item Number: CM0026725



Introduction

Rabbit polyclonal antibody targeting human NY-ESO-1 (CTAG1), validated for WB, IHC-P, ELISA. Ideal for cancer/testis antigen research, as NY-ESO-1 is expressed in various cancers and testis. Produced against amino acids 1-100.

[Learn More](#)

Parameter	Value
Product Name	NY-ESO-1 Rabbit pAb
Remarks/Alias	CTAG; ESO1; CT6.1; CTAG1; LAGE-2; LAGE2B; NY-ESO-1
Species	Human
GeneID (Human)	1485
GeneID	1485
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human NY-ESO-1 (NP_001318.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P78358
Protein Weight	18kDa
Shipping	Ice bag