

Anti-Ny-Eso-1 Rabbit Monoclonal Antibody For Wb, Elisa - P78358

Item Number: CM0002683



Introduction

Rabbit monoclonal antibody targeting human cancer/testis antigen NY-ESO-1 (CTAG1). Validated for WB and ELISA applications, detects an 18kDa band, suitable for oncology and immunodiagnostic research.

[Learn More](#)

Parameter	Value
Product Name	NY-ESO-1 Rabbit mAb
Remarks/Alias	CTAG; ESO1; CT6.1; CTAG1; LAGE-2; LAGE2B; NY-ESO-1
Species	Human
GeneID (Human)	1485
GeneID	1485
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human NY-ESO-1 (P78358)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P78358
Protein Weight	18kDa
Shipping	Ice bag